**KIT ANALISI**

**ANALYSIS KIT**

**36**

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## MEASUREMENT KITS

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**THE USE OF HAZARDOUS SUBSTANCES FOR THE CHEMICAL ANALYSIS OF WATER**

Water is a fundamental resource for life and the industry. The chemical analysis are required to establish its purity before it can be utilized for personal hygiene, alimentary uses and industrial purposes, because considering its vital and final use should not contain anything hazardous for all the living species and for all the industrials' machinery. Even the waisted water should be appropriately analyzed and purified to protect the environment and to avoid the contamination of the natural resources starting to become lacking respect the overall needs. To analyze the majority of chemical substances several methods are available such as potentiometry, polarography, different chromatography, gravimetry, volumetry, several photometry etc.. Generally the fastest methods, also more reliable and low cost are suitable for analysis in situ, but apply chemical reagents whereas the methods that implies expensive machinery, usually not available as a transportable unit, do not apply chemical reagents or their use is rather limited. Often the results are required to be obtained immediately demanding a portable and limited equipment to allow all the analysis in situ. Reason that has contributed to make the easier and cheaper methods to become more utilized due to their overall reliability and fast results achieved for example the titration and the spectrophotometric analysis that often can be performed by employees not specific qualified and trained. These latest mentioned analysis are based on the use of a discrete amount of chemical reagents, that in some instances, are classified among the substances hazardous for health. The modern legislation requires that every chemical product utilized in the industry environment needs a technical risk report clarifying its properties, pointing out any possible correlated risks and giving some advice on the security procedures to adopt, also all the hazardous substances used should be clearly identified. These overmentioned technical risk reports, without any doubt useful to everyone, can sometime generate in an individual without any chemical skills unjustified fears that can affect the use of chemical reagents indispensable for certain analysis, even if their use is substantially proved to be without any health risk for the user. An example: a bottle containing a reagent with 40% of Ethyl Alcohol in a industry environment needs to be clearly labelled as "high flammable" and is perceived as dangerous, whereas, in an house, a 1.000 ml bottle of any sort with the same Alcohol concentration or more does not arise any worries (the actual Law for perfumes or spirits does not require a "high flammable" label). The legislator has to resolve the following dilemma. He has, to one side, to safeguard the employees health and safety minimizing their direct or indirect work exposure. At the other side a wrong or an extreme interpretation of certain dispositions could lead to

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to create unreal problems, mostly if the final readers has no adequate technical skills, more evident in the chemical environment. Enlightening is the research conducted by Roger Bate, a economist from Cambridge, several years ago that involved interviewing common people. He asked the following: "Are you favorable to the abolish of the Dihydrogen monoxide?" 76% of all the interviewed answered affirmatively. Their cultural level was not clearly stated, but is a fact that to a non chemist the expression "Dihydrogen monoxide" is not immediately linked to a substance composed of two Hydrogens and a single Oxygen, that is H<sub>2</sub>O - water. This research clearly demonstrate that the sensibilization of the common peolple towards environmental problems, laudable conducted by the local authority and by the mass media, can sometimes unintentionally generate real cases and seriuos problems of "ambiental hysteria" that sometimes has the potential risk of contaminating even the thoughts of the people called to check and control the working environment and the same substances required to complete the working task. A serious problem, that can rise from a wrong law interpretation, is that many people consider the use of any kind of dangerous substance illegal or must be abolished at any cost. This is incorrect and not true. If the production of a certain number of substances potentially hazardous is permitted, is simply because considering the actual technical level our society is unable to do without them. Simply by considering the amount of poisons involved to defeat insects or kill bugs, or to all the amount of medicine used, or to all the preservatives that prevent the entire harvesting that will potentially decay without them.

Obviously the use of this kind of substances is legal, only if the preventive measurements are implemented to minimise their connected risks. At this point another problem emerge in the industry where this dangerous products are manufactured and manipulated, the employees are continuously exposed to hight product concentrations, here the substances need to be loaded and unloaded from the industrial machinary and from the silos or into storage containers, at the end into their final packeging ready to be sold or to be used, whereas in the factory that only use the final product the employees exposure risks are quite limited. In some extreme cases dangerous reagents enter the final products at a diluited concentration, as in the case of the reactants used for the analysis of water in solutions all ready for use. The exposure risk for employees that irregularly handle the final product is irrelevant, as in the overmentioned reactants, especially if the final reactants are contained in little bottles with water drip end that can be directly mixed with the water samples that needs to be analyzed. The products classfied as hazardous according with the D.L.03/02/97 N. 52 are approximately 3.500, but the cases considered are more numerous due to the level of risk that for several products is directly proportional to their effective final concentration. Is quite obvious that the products attached risk reports can not contemplate several levels of

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hazard resulted from different improper uses, but they must provide enough information that will allow any employee to work in a safe environment and in any critical situation. Its quite comprehensible that in the risk report regarding Sodium Chloride, known as common salt, textually reports: "6. Procedures in case of any accidental spillage: Wear an individual respiratory, safety goggles, rubber boots and resistant rubber gloves. Sweep.....etc." The protective use of gloves and safety goggles is undoubtedly adequate for an employee that works in a salt mine, but appears exaggerated for a person that sprinkles a pinch of salt over a salad plate. Substances that we habitually ingest, like salt, or that are vital for our survival, like water, is a matter of fact that if they are ingested in an excessive quantity they can be dangerous for our health or even sometimes become lethal. Without wanting to interfere with the validity of the identification and classification of hazardous substances, its necessary to distinguish as well as, or even between SITUATIONS MORE OR LESS DANGEROUS and bend the safety procedures in fuction of them. If we want to use the electric light we must protect ourselves from any eventual electric discharges; instead if we wish to use a light of a candle we must prevent the potential risk of fire, in any case the employee dealing with the production of electric energy, or with the maintenance of hight voltage lines should wear an adequate personal protective clothing and equipment respect a person that reads using an electric bulb. The law maker or the public or private employer must take, without any doubt, all the necessary precautions to minimise any possible risks for the employee whilst dealing with hazardous substances. Ideally, in order to achieve this, it will be ideal to eliminate all the potential dangers from any working activities, but unfortunately this is not always possible, in the same way its impossible to prevent the risk of fire for either a car driver or an air pilot but even for a bicycle rider or for a person that simply walks in the street. The risk reports are created with the aim of providing enought information to the competent individual in order to enable him to assess the potential risks and work accordingly to prevent any harms to the workers, but these reports handled by persons without a particular competence for making these kind of decisions can result counter-productive interfering with the normal course of the work. For this the directive 626/94 - 242/96 "Workers Health and Safe" sets precise competences for each professional figure. The employer, the manager, the security responsible are all involved in the risks assesment, in the choice of the equipments and with the workers training and education. Regarding the water analysis its necessary to distinguish between two distinctive fundamental situations. The first where the analysis are performed into a set-up laboratory where the indispensable reagents are created directly by the workers by weighing, mixing and diluting the primary materials (pure chemical products for analysis). In this instance the workers must have

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an adequate education ( i.e. A Chemistry degree or diploma or being an analyst ), that includes how to manipulate hazardous chemical products as part of their training. The risks are the normal ones for a chemical laboratory. The risk reports can be handed to trained workers without any further comment. In the second case the analysis are performed by untrained workers by means of commercial analytical kits, that contain reagents ready for use. In this instance the workers are not required to handle hazardous substances, but they only need to add directly to the water samples a discrete amount of drops of the reagents all ready for use ( or a microspoon of powder ), otherwise inject a few milliliters of the water sample into a vial already coated with the active reagent. In this last case the potential risk is negligible but the safe environment should be the same assessed by a person with an adequate competence, that has also the duty to instruct the workers providing them with further explanations in order to avoid a misreading that can cause several problem. This audit trail becomes more essential when often several regulations discrepancy are encountered or when to decide what will be the best analysis to perform in a certain situation. Considering the following example: the recommended method to analyze Chlorine contained in water is by using a colorimeter with O-Tolidine. Being O-Tolidine a carcinogenic product, the chemical industry has created an alternative method that implies the reagent DPD, this has the advantage to be able to measure the Chlorine mixed or total (whereas with O-Tolidine its only possible to determine the total).

The reagent DPD is in turn toxic and more complicated, the worker or the person in charge should be able to make an informed decision based, as well as, on the result desired. A further example is: the method specified for the analysis of COD requires the use of a 50 ml water sample with the addition of bichromate of potassium, sulphuric acid concentrated and silver sulphate and after boiling in an appropriate glass apparatus, a titration is followed with ammonium sulphate and iron. The overmentioned analysis should be performed into a laboratory by workers with chemical skills. A more modern method, but not officially recognized, requires the use of samples in the range of 0.5 to 2 ml, being 25 to 100 times less. The reagents quantities, proportionally more reduced, are placed into pre-weighted vials ready for use, in which the worker only need to add the water sample. The measurement in this instance it is not done by means of titration but through colorimetric analysis applying a spectrophotometer. This method doesn't require any reagents manipulation and for this reason its application can be relied on untrained personnel, but after having an adequate level of training.

In case of any doubt, the assessment of any potential or real hazardous situation should be performed by an individual with appropriate skills that is able to correctly interpret all the technical risk reports without superficiality and even without fanaticism. To clarify this concept consider the following example: one of

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the reagents most utilized in the analysis of the water is the hardness indicator, that often contains triethanolamine, considered an "dangerous ingredient". A chemical graduate knows that is a product utilized for example in the cosmetic industry and is able to evaluate the potential risk for an operator dealing with this "dangerous ingredient" in case of any drops accidentally come into contact with his skin.

Whereas an operator without any chemical education probably is unaware that a normal detergent or a disincrustant that he normally use, without any worries, at home is potentially more dangerous leading him to refuse to adopt the kit to determine the water hardness or making him to adopt excessive protective measurements that can obstruct his normal duty only because frightened by the relative risk technical report. Concluding we can establish the following:

- The preparation of the reagents by means of manipulating the relative primary materials must be done exclusively by qualified operators working in a chemical laboratory
- The reagents contained in the ready for use vials part of the analytical kits in commerce can be used even by untrained personnel
- How uses the analytical kits must follow the instructions, along with the hygienic and working norms prescribed by Law
- The risk assessment technical reports regarding the reagents should be meticulously valued by qualified personnel and where necessary explained and expounded towards the untrained personnel to avoid any inconsiderate use or preventing the rise of groundless fears
- The commercial analytical kits have been in use since several years without having created any health and safety issues
- The contact between the operator and the reagents is highly unlikely, except in cases of serious negligence (i.e. analysis performed instead of in its appropriate container, in a common glass used later on to drink without washing it previously) or due to the breaking of a vial or vessel containing a reagent (in this case the amount of reagent spilled and the possible skin contact are however well under the hazard limits). Substantially the risk reports are useful especially in these kind of situation providing advice on how to handle an unexpected event
- The water analysis done with the commercial analytical kit do not present any risk for the operator that follows diligently the instructions provided and the rules dictated by the good sense before applying the actual hygiene at work legislation.

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## ACIDITY KIT

### OBSERVATIONS:

- The classical indicator (methyl orange)from which the analysis gather is name has been replaced by the "indicator m", that has a more accurate toning, because it allows to identify the final titration point more easily.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ Methyl orange acidity
- METHOD: \_\_\_\_\_ Sodium hydrate titration
- SENSIBILITY: \_\_\_\_\_ 0.5 French degrees (on a 10 ml sample)
- DETERMINATIONS NUMBER: \_\_\_\_\_ 625 French degrees

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## ALKALINITY KIT

### OBSERVATIONS:

The classical indicator (methyl orange) from which the analysis gather is name has been replaced by the "indicator m", that has a more accurate toning, because it allows to identify the final titration point more easily.

- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ Methyl orange alkalinity or phenolphthalein in °f ( created for low alkalinity waters)
- METHOD: \_\_\_\_\_ Hydrochloric acid titration
- SENSIBILITY: \_\_\_\_\_ 0.5 French degrees (on a 10 ml sample)
- DETERMINATION NUMBER: \_\_\_\_\_ 625 French degrees

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## ALKALINITY BOILERS KIT

### OBSERVATIONS:

- The classical indicator (methyl orange) from which the analysis gather is name has been replaced by the "indicator m", that has a more accurate toning, because it allows to identify the final titration point more easily.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ Alkalinity using methyl orange and phenolphthalein in epm
- METHOD: \_\_\_\_\_ Hydrochloric Acid titration
- SENSIBILITY: \_\_\_\_\_ 0.5 epm (10 ml sample)
- DETERMINATIONS NUMBER: \_\_\_\_\_ 625 epm

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
### M ALKALINITY KIT

**OBSERVATIONS:**

- The classical indicator (methyl orange)from which the analysis gather is name has been replaced by the "indicator m", that has a more accurate toning, because it allows to identify the final titration point more easily.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ Methyl orange alkalinity in °f ( created for low methyl orange alkalinity waters, such as untreated or softened waters. The german habit is to identify the alkalinity m with the "carbonatic hardness",meaning temporary hardness. This is true only in absence of sodium bicarbonate. )
- METHOD: \_\_\_\_\_ Hydrochloric acid titration
- SENSIBILITY: \_\_\_\_\_ 0.5 French degrees (on a 10 ml sample)
- DETERMINATION NUMBER: \_\_\_\_\_ 625 French degrees

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
## P ALCALINITY KIT

### OBSERVATIONS:

- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ Phenolphthalein alkalinity in epm (ideal for high phenolphthalein alkalinity waters)
- METHOD: \_\_\_\_\_ Hydrochloric Acid titration
- SENSIBILITY: \_\_\_\_\_ 0.5 epm (10 ml sample)
- DETERMINATIONS NUMBER: \_\_\_\_\_ 625 epm

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
### TOTAL ALKALINITY KIT

**OBSERVATIONS:**

- The classical indicator (methyl orange)from which the analysis gather is name has been replaced by the "indicator m", that has a more accurate toning, because it allows to identify the final titration point more easily.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ Methyl orange alkalinity in epm (created for high total alkalinity waters, such as waters of generators of vapor)
- METHOD: \_\_\_\_\_ Hydrochloric Acid titration
- SENSIBILITY: \_\_\_\_\_ 0,5 epm
- DETERMINATION NUMBER: \_\_\_\_\_ 625 epm

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**CERTIFICATI:**

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**AMMONIA KIT****OBSERVATIONS:**

- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ COLORIMETRICAL

- METHOD: \_\_\_\_\_ Nessler

- SENSIBILITY: \_\_\_\_\_ 0,25 ppm NH<sub>3</sub>

- DETERMINATION NUMBER: \_\_\_\_\_ 200

**KIT CONTENT:**

- 1 bottle 25cc reagent Ammonia A
- 1 bottle 25cc reagent Ammonia B
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- instructions

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
### DIOXIDE CARBON KIT

**OBSERVATIONS:**

- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ Phenolphthalein acidity
- METHOD: \_\_\_\_\_ Sodium hydrate titration
- SENSIBILITY: \_\_\_\_\_ 0,5 French degrees (on a sample of 10 ml)
- DETERMINATION NUMBER: \_\_\_\_\_ 625 French degrees

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## STRONG ANIONS KIT

### OBSERVATIONS:

- Using the anions strong activator, these are transformed in strong acids and titrated as acids. For this reason the kit **must be combined with the IPT Acidity kit m**

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ Transformation in strong acids followed by methyl orange acid titration
- METHOD: \_\_\_\_\_ Titration with sodium hydroxid
- SENSIBILITY: \_\_\_\_\_ 0,5 French degrees (on a 10 ml sample)
- DETERMINATION NUMBER: \_\_\_\_\_ approx 30

### KIT CONTENT:

- 50 ml Activator strong anions (into 2 containers 25ml each)
- 1 bottle 25cc with water drip, empty
- 1 graduated syringe (without piston)
- 1 spoon N.4
- instructions

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
## CALCIUM KIT

### OBSERVATIONS:

- The indicator toning is clearly visible ( in opposition to the old indicator toning, that required a trained eye)
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ Volumetric ( created to complete the hardness analysis, in cases where the hardness percentage related to Calcium or Magnesium is required )
- METHOD: \_\_\_\_\_ Calcium titration with EDTA
- SENSIBILITY: \_\_\_\_\_ 0,5 French degrees ( on a 10 ml sample
- DETERMINATION NUMBER: \_\_\_\_\_ 625 French degrees

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## TOTAL CHLORINE

### OBSERVATIONS:

- The O-tolidine concentration in the reagent total Chlorine is lower than the limit established by the directive 67/548 CEE to be classified as dangerous
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric, created for chlorate waters with few organic substances ( for example drinkable waters)
- METHOD: \_\_\_\_\_ O-tolidine
- SENSIBILITY: \_\_\_\_\_ 0,1 ppm Cl<sub>2</sub>
- DETERMINATION NUMBER: \_\_\_\_\_ 75

### KIT CONTENT:

- 1 bottle 25cc reagent total Chlorine
- 1 test-tube from 10cc with lid
- 1 chromatic scale 0,1-1,5 ppm
- instructions

### ORDERING CODE:

- HC0011

### REPLACEMENTS AND ACCESSORIES:

- HC0122: reagent total Chlorine, 6 bottles 25cc per pack
- HC0116: test-tube 10cc, 12 tubes per pack

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### DPD CHLORINE KIT

#### OBSERVATIONS:

- This kit allows to determine both the free chlorine and the chlorine total. The difference between the two values corresponds to the chlorine combined
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

#### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric (created for waters with organic substances, i.e. the water of the swimming pool)-
- METHOD: \_\_\_\_\_ n,n-diethyl-phenyldiamine
- SENSIBILITY: \_\_\_\_\_ 0,1 ppm Cl<sub>2</sub>
- DETERMINATION NUMBER: \_\_\_\_\_ 150

#### KIT CONTENT:

- 1 bottle 15cc reagent Chlorine DPD A
- 1 bottle 15cc reagent Chlorine DPD B
- 1 bottle 15cc reagent Chlorine DPD C
- 1 test-tube 10cc with lid
- 1 chromatic scale 0,1-1,5 ppm
- instructions

#### ORDERING CODE:

- HC0012

#### REPLACEMENTS AND ACCESSORIES:

- HC0123: reagent Chlorine DPD A, 6 bottles 15cc per pack
- HC00124: reagent Chlorine DPD B, 6 bottles 15cc per pack
- HC00125: reagent Chlorine DPD C, 6 bottles 15cc per pack
- HC0116: test-tube 10cc, 12 tubes per pack

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## TITRANT CHLORIDES KIT

### OBSERVATIONS:

- Considering the "threshold" of the methods of determination of the chlorides with titration, we advise the colorimetric method (cod.HC0014) for low concentrations
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ VOLUMETRICAL
- METHOD: \_\_\_\_\_ Titration with nitrate mercuric
- SENSIBILITY: \_\_\_\_\_ 0,5 French degrees (on a 10 ml sample)
- DETERMINATION NUMBER: \_\_\_\_\_ 625 French degrees

### KIT CONTENT:

- 1 bottle 25cc titrant chlorides
- 1 bottle 15cc reagent chlorides A
- 1 bottle 15cc reagent chlorides B
- 1 test-tube from 20cc
- instructions

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### COLORIMETRIC CHLORIDES KIT

**OBSERVATIONS:**

- In the range between 0 and 80 ppm Cl the colorimetric method is the most titration sensible.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric
- METHOD: \_\_\_\_\_ Sulpho-cyanide ferric
- SENSIBILITY: \_\_\_\_\_ 1 ppm Cl
- DETERMINATION NUMBER: \_\_\_\_\_ 80

**KIT CONTENT:**

- 1 bottle 25cc reagent Chlorides A col.
- 1 bottle 25cc reagent Chlorides B col.
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- instructions

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## COLORIMETRIC CHROMED KIT

### OBSERVATIONS:

- The method detects the Cr<sup>+6</sup>. To detect the Cr<sup>+3</sup> is necessary to detect the chromium total using the laboratory photometric method (see Hydrocheck photometer manual). Il Cr<sup>+3</sup> is obtained from the difference between chromium total and Cr<sup>+6</sup>.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric
- METHOD: \_\_\_\_\_ Diphenilcarbazide
- SENSIBILITY: \_\_\_\_\_ 0,05 ppm Cr
- DETERMINATION NUMBER: \_\_\_\_\_ 50

### KIT CONTENT:

- 1 bottle 25cc reagent Chromate A
- 1 bottle of powder reagent Chromate B
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- 1 small spatle
- instructions

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### DEMINERALIZATION KIT

**OBSERVATIONS:**

• The aim of the kit is to rinse the ion, cation and anion exchanges filters that make a demineralization system. After the regeneration the cationic effluent has an high acid concentration. The rinses should be protracted until the pH indicator is at a constant value between pH 0-5. In the same way, after the regeneration the anionic effluent has an elevated concentration of alkali. The rinses should be protracted until the pH indicator is at a constant value between pH 1-11.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric (ideal for the control of demineralization systems)
- METHOD: \_\_\_\_\_ determination of the pH in the acid and alkaline fields.
- SENSIBILITY: \_\_\_\_\_ 0.5 pH units

**KIT CONTENT:**

- 1 bottle 25cc reagent pH 1-11
- 1 bottle 15cc reagent pH 0-5
- 1 test-tube 10cc with lid
- instructions

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## HARDNESS KIT

### OBSERVATIONS:

- Strictly together with the Calcium and Magnesium are titrated even the Strontium and Barium, nevertheless the concentration of the two latest ones is so lower that can be neglected.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ volumetric
- METHOD: \_\_\_\_\_ Metals alkaline-earthly titration ( above all Calcium and Magnesium) with EDTA
- SENSIBILITY: \_\_\_\_\_ 0,5 French degrees ( on a 10 ml sample)
- DETERMINATION NUMBER: \_\_\_\_\_ 625 French degrees

### KIT CONTENT:

- 1 bottle 25cc hardness titration
- 1 bottle 15cc hardness indicator
- 1 test-tube from 20cc
- instructions

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## HARDNESS BOILERS KIT

### OBSERVATIONS:

• In cases where the thermic circuit is proceeding, that includes a systematic control of the water hardness that feeds the vapour generator, and the maintenance of a pad against the water hardness (for example phosphate) in the boiler, is not necessary to determine the hardness of the water contained in the vapour generator. Besides substances present in the water from the boiler can interfere this kind of analysis, by inhibiting the toning of the hardness indicator and by consequently simulating the presence of hardness, that in reality there is not. Wanting to analyze the water hardness, nevertheless, is recommended to add the hardness reagent A, that in several cases, it is able to eliminate this interference. The hardness reagent A is not usable for waters that have a very low conductivity or with a pH lower than 9.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ volumetric
- METHOD: \_\_\_\_\_ Hardness titration with inhibition of some interferences.
- SENSIBILITY: \_\_\_\_\_ 0,5 French degrees ( on a 10 ml sample)
- DETERMINATION NUMBER: \_\_\_\_\_ 625 French degrees

### KIT CONTENT:

- 1 bottle 25cc hardness titration
- 1 bottle 15cc hardness indicator
- 1 bottle 15cc hardness reagent A
- 1 test-tube from 20cc
- instructions

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## RESIDUAL HARDNESS KIT

### OBSERVATIONS:

- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric  
(ideated to reveal any traces of hardness)

- METHOD: \_\_\_\_\_ Identify traces of hardness by mean of an indicator that assumes different colorations in the range of 0 and 3 ppm.

- SENSIBILITY: \_\_\_\_\_ 0,25 ppm CaCO<sub>3</sub>

- DETERMINATION NUMBER: \_\_\_\_\_ 250

### KIT CONTENT:

- 1 bottle 25cc reagent Hardness A
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- instructions

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## IRON KIT

### OBSERVATIONS:

- The chromatic scale 0,25-15 ppm is utilized with a 5 ml sample. To obtain a higher sensibility the chromatic scale 0,05 ppm is used on a 20 ml sample ( increasing in proportion the reagent drops)
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetrica
- METHOD: \_\_\_\_\_ 1,10-fenantrolina
- SENSIBILITY: \_\_\_\_\_ 0,05 ppm Fe
- DETERMINATION NUMBER: \_\_\_\_ 100 (with scale 0,25-15 ppm)

### KIT CONTENT:

- 1 bottle 25cc reagent Iron A
- 1 test-tube from 10cc with lid
- 1 test-tube from 20cc with lid
- 1 chromatic scale 0,25-15 ppm
- 1 chromatic scale 0,05-0,2 ppm
- instructions

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
### 1 - 20 ppm PHOSPHATE KIT

#### OBSERVATIONS:

- The method identify the Phosphorus in the form of orthophosphate.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

#### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric
- METHOD: \_\_\_\_\_ molybdate
- SENSIBILITY: \_\_\_\_\_ 1 ppm PO<sub>4</sub>
- DETERMINATION NUMBER: \_\_\_\_\_ 100

#### KIT CONTENT:

- 1 bottle 25cc reagent Phosphate A
- 1 bottle 25cc reagent Phosphate B
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- instructions

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
**CERTIFICATI:**

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**OBSERVATIONS:**

- The method reveals the presence of Phosphorus in the form of orthophosphate. Generally the IPT phosphate kit 1-20 are preferred because their chromatic scale has a better color variation.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric
- METHOD: \_\_\_\_\_ molybdate vanadate to determine elevated quantity of orthophosphate
- SENSIBILITY: \_\_\_\_\_ 5 ppm PO<sub>4</sub>
- DETERMINATION NUMBER: \_\_\_\_\_ 50

**KIT CONTENT:**

- 1 bottle 15cc reagent Phosphate A
- 1 bottle 25cc reagent Phosphate B
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- instructions

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## HYDRAZINE KIT

### OBSERVATIONS:

- Using the classical method even the samples without Hydrazine have a slightly yellow coloration. Whereas with our modified method samples without Hydrazine remain colorless. Therefore the light yellow coloration identified using our method corresponds to an effective presence of Hydrazine in traces.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetrica
- METHOD: \_\_\_\_\_ p-dimethylaminobenzaldehyde
- SENSIBILITY: \_\_\_\_\_ 0,1 ppm N<sub>2</sub>H<sub>4</sub>
- DETERMINATION NUMBER: \_\_\_\_\_ 50

### KIT CONTENT:

- 2 bottles 25cc reagent Hydrazine A
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- instructions

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**CERTIFICATI:**

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**OBSERVATIONS:**

- In the water depuration arena the terms "nitrogen nitrous" (that correspond to nitrates), and "nitrogen ammoniacal"(that refers to ammonia)are often quoted, nevertheless usually the concentrations are often referred as NO<sub>3</sub> for the nitrates, as NO<sub>2</sub> for i nitriti and as NH<sub>3</sub> for ammonia.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ COLORIMETRICAL
- METHOD: \_\_\_\_\_ Acid 5-hydroxyallicic
- SENSIBILITY: \_\_\_\_\_ 10 ppm NO<sub>3</sub>
- DETERMINATION NUMBER: \_\_\_\_\_ 50

**KIT CONTENT:**

- 1 bottle powder reagent Nitrates A
- 1 test-tube from 10cc with lid
- 1 cromatic scale
- 1 spoon
- instructions

**ORDERING CODE:**

- HC0024

**REPLACEMENTS AND ACCESSORIES:**

- HC0140: reagent A, 6 bottles per pack
- HC0116: test-tube 10cc, 12 tubes per pack

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## NITRITE KIT

### OBSERVATIONS:

• In the water depuration arena the terms "nitrogen nitrous" (that correspond to nitrates), and "nitrogen ammoniacal" (that refers to ammonia) are often quoted, nevertheless usually the concentrations are often referred as NO<sub>3</sub> for the nitrates, as NO<sub>2</sub> for i nitriti and as NH<sub>3</sub> for ammonia.

• The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetrica
- METHOD: \_\_\_\_\_ Sulfanilic Acid
- SENSIBILITY: \_\_\_\_\_ 0,05 ppmNO<sub>2</sub>
- DETERMINATION NUMBER: \_\_\_\_\_ 50

### KIT CONTENT:

- 1 bottle 25cc reagent Nitriti A
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- instructions

### ORDERING CODE:

- HC0025

### REPLACEMENTS AND ACCESSORIES:

- HC0141: reagent A, 6 bottles per pack
- HC0116: test-tube 10cc, 12 tubes per pack

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**CERTIFICATI:**

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**OBSERVATIONS:**

- Normally it is enough to analyze waters using the universal pH kit, with a range from 1 to 11. The kit 0-5 is only adopted when the analysis require a higher level of precision.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ colorimetrica
- METHOD: \_\_\_\_\_ Mixed indicator
- SENSIBILITY: \_\_\_\_\_ 1 pH unit
- DETERMINATION NUMBER: \_\_\_\_\_ 200

**KIT CONTENT:**

- 1 bottle 25cc indicator pH 0-5
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- instructions

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### PH KIT 6.8 - 8.4

#### OBSERVATIONS:

- The kit offers an high level of precision in the field of great interest from a level lightly acid to a level lightly alkaline (6.8-8.4), thanks to the mixed indicator, that in this range offers a various color changes ( from yellow to violet, passing through the orange and red) therefore allows to appreciate even very small differences of pH with a certain precision.
- A more accurate result can be achieved using an electric pH-meter of good quality

#### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric (to analyze the pH in swimming-pools and in thermal circuits)
- METHOD: \_\_\_\_\_ mixed indicator
- SENSIBILITY: \_\_\_\_\_ 0,2 pH units
- DETERMINATION NUMBER: \_\_\_\_\_ 200

#### KIT CONTENT:

- 1 bottle 25cc indicator pH 6.8-8.4
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- instructions

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**CERTIFICATI:**

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**OBSERVATIONS:**

•The selection indicator is ideal for all the cases where a greatest precision is not required (obtainable instrumentally with a good pH-meter).In the most important field, from 5 to 10, is possible a good color differentiation that varies from 0,5 to 0,5 pH units, whereas in the field lower than 5 the pH differences varies from unit to unit and they involve a minimal color differentiation in the red, difficult to interpret even for a well trained eye. For this reason has been created the pH 0-5 kit (cod. 1020), that offers a better color differentiation in this range. Besides in the interesting field from slightly acid to slightly alkaline (6.8-8.4) the pH 6.8-8.4 kit has been created (cod.1021) with a mixed indicator that, in this field, offers a great differences of color (from yellow to violet, passing through orange and red) and therefore allows to appreciate ,with a certain degree of precision,even small differences of pH.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric
- METHOD: \_\_\_\_\_ mixed indicator for an "universal" pH determination
- SENSIBILITY: \_\_\_\_\_ 0,5 pH units
- DETERMINATION NUMBER: \_\_\_\_\_ 200

**KIT CONTENT:**

- 1 bottle 25cc indicator pH 0-11
- 1 test-tube from 10cc with lid
- 1 chromatic scale
- instructions

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## SWIMMING POOLS KIT

### OBSERVATIONS:

• Considering the fact that swimming pools water contain organic substances that can combine with Chlorine forming chloramines, is recommended to perform the analysis of free Chlorine, of total Chlorine and combined Chlorine, using the DPD method.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetric
- METHOD: \_\_\_\_\_ Measurement of the pH using the mixed indicator and of the Chlorine using the DPD method
- SENSIBILITY: \_\_\_\_\_ 0,2 pH units/0,1 ppm Cl<sub>2</sub>
- DETERMINATION NUMBER: \_\_\_\_\_ 120/150

### KIT CONTENT:

- 1 bottle 15cc reagent Chlorine DPD A
- 1 bottle 15cc reagent Chlorine DPD B
- 1 bottle 15cc reagent Chlorine DPD C
- 1 bottle 15cc reagent pH 6.5-8.5
- 1 test-tube from 10cc with lid
- 1 syringe 5cc graduated
- 1 chromatic scale
- instructions

### ORDERING CODE:

- HC0029

### REPLACEMENTS AND ACCESSORIES:

- HC0123: reagent A, 6 bottles 15cc per pack
- HC0124: reagent B, 6 bottles 15cc per pack
- HC0125: reagent C, 6 bottles 15cc per pack
- HC0091: reagent pH, 6 bottles 15cc per pack
- HC0116: test-tube 10cc, 12 tubes per pack

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## SCALING HARDNESS KIT

### OBSERVATIONS:

• This kit is usually used to demonstrate the difference between water conditioned and non conditioned. In water without hardness there is no precipitate. Whereas in water with hardness there is the formation of precipitate that is as more consistent as the level of hardness, even if it doesn't provide an analytical measurement of the hardness quantity.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ Demonstration of the presence of hardness through formation of a precipitate.

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## SILICA KIT

### OBSERVATIONS:

- This method is very sensible and allows to determine concentration even lower than those correctly indicated. Result achieved by simply observing the color formed on the top layer of 5 - 10 cm, contained into a transparent, white base container adjusting the number of reagents drops to the quantity of sample used.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ colorimetrica
- METHOD: \_\_\_\_\_ molybdate
- SENSIBILITY: \_\_\_\_\_ 0,25 ppm SiO<sub>2</sub>
- DETERMINATION NUMBER: \_\_\_\_\_ 50

### KIT CONTENT:

- 1 bottle 15cc Silica reagent A
- 1 bottle 15cc Silica reagent B
- 1 bottle 15cc Silica reagent C
- 1 test-tube from 10cc with lid
- 1 chromatic scale 0,25-8 ppm
- instructions

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### SULFITE KIT 2 ppm

#### OBSERVATIONS:

- In the IPT sulfite kit 2 ppm (cod.HC0032), every titration drop employed on a 5 ml sample corresponds to 2 ppm SO<sub>3</sub>. In the case where is necessary to identify relatively high concentrations of sulfite, is possible to use the sulfite kit 10 ppm (cod. HC0033), inn which every titration drop employed on a 5 ml sample corresponds to 10 ppm SO<sub>3</sub>.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

#### CHARACTERISTICS:

- ANALYSIS TYPE: \_\_\_\_\_ volumetrical
- METHOD: \_\_\_\_\_ potassium iodided-iodate titration for an "universal" determination
- SENSIBILITY: \_\_\_\_\_ 1 ppm SO<sub>3</sub> (on a 10 ml sample)
- DETERMINATION NUMBER: \_\_\_\_\_ 1.250 ppm SO<sub>3</sub>

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### SULFITE KIT 10 ppm

#### OBSERVATIONS:

- In the IPT sulfite kit 2 ppm (cod. HC0032), every titration drop employed on a 5 ml sample corresponds to 2 ppm SO<sub>3</sub>. In the case where is necessary to identify relatively high concentrations of sulfite, is possible to use the sulfite kit 10 ppm (cod. HC0033), in which every titration drop employed on a 5 ml sample corresponds to 10 ppm SO<sub>3</sub>.
- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

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**CERTIFICATI:**

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**OBSERVATIONS:**

- The IPT kit are ideal in cases where is necessary to perform only a few analysis with a discrete precision.

**CHARACTERISTICS:**

- ANALYSIS TYPE: \_\_\_\_\_ colorimetrical-
- METHOD: \_\_\_\_\_ Formaldossime
- SENSIBILITY: \_\_\_\_\_ 0,1 ppm Mn
- DETERMINATION NUMBER: \_\_\_\_\_ 70

**KIT CONTENT:**

- 1 bottle 15cc reagent Manganese A
- 1 bottle 15cc reagent Manganese B
- 1 bottle 15cc reagent ManganeseC
- 1 test-tube from 10cc with lid
- 1 cromatic scale
- instructions

**ORDERING CODE:**

- HC0034

**REPLACEMENTS AND ACCESSORIES:**

- HC0137: reagent A, 6 bottles 15cc per pack
- HC0138: reagent B, 6 bottles 15cc per pack
- HC0139: reagent C, 6 bottles 15cc per pack
- HC0116: test-tube 10cc, 12 tubes per pack

**CERTIFICATES:**

- Health and Safety risk assesement